

## RoweEpiflow



- According to new standard DIN EN ISO 80369 part 6: „Connectors for neuraxial applications“
- 20% smaller than Luer connections
- Prevents risk of misconnections
- Visual drug flow control
- For epidural and neuraxial anaesthesia
- Male connector with rotating lock nut (yellow color coding)
- Without venting membrane

Art.-No.	Material (body)	Filling volume (ml)	Connection	Filter	PU ((pcs./box) / pal.)
A-6750	PC	0.4	female / male	0.2 µm PET	1,000 / 16,000
A-6751	PC	0.4	female / male	0.2 µm Nylon *	1,000 / 16,000
A-6752	PC	0.4	female / male	0.2 µm PES	1,000 / 16,000

\* positively charged

RoweMed AG - Medical 4 Life  
 www.rowemed.de • info@rowemed.de  
 Juri-Gagarin-Ring 4 • 19370 Parchim • Germany  
 Tel.: +49 (0) 3871 451280 • Fax : +49 (0) 3871 451282

